

EXCLUSIVE TO SLS

NEW PIPETGIRL

Let's fight breast cancer together!



A £10 donation per PIPETGIRL goes to breast cancer research institutes



Our commitment to support breast cancer research

Since 1965, INTEGRA Biosciences has been dedicated to developing innovative liquid handling tools for the Life Science industry.

It is our passion to provide scientists with reliable instruments of the highest quality and precision, so they can focus all their efforts on research.

On our 50 year anniversary, we would like to support cancer research by donating £10 for every purchase of a PIPETGIRL pipette controller to research institutes dedicated to finding a cure for breast cancer.



**SCIENTIFIC
LABORATORY
SUPPLIES**

INTEGRA

www.integra-biosciences.com

NEW PIPETGIRL

PIPETGIRL offers the same features as INTEGRA's PIPETBOY acu 2: an increased pipetting speed, a significant weight reduction and unmatched pipetting precision.



- Lighter, with only 195 grams and an optimized ergonomic balance.
- Faster, the new Turbo Mode provides the unit with 20% extra pipetting speed.
- Longer operating time due to the latest lithium ion battery technology.

Instrument	Part No.
PIPETGIRL	I155 021
Complete with battery Lithium Ion, mains adapter, spare sterile filter 0.45 µm and wall mount	

Accessories	Part No.	
Stand	for PIPETGIRL and VACUBOY	I155 065
Wall mount		I155 521
Mains adapter (100-240 VAC, 50/60 Hz)	UK version	I153 214

Consumables	Part No.	
Filter	0.45 µm sterile	I153 015
	0.45 µm unsterile	I153 016
	0.2 µm sterile	I156 608
	0.2 µm unsterile	I156 607
Pipette mount	silicone, for holding the pipette in the sterile module	I151 020
Filter rubber	silicone, for holding the filter in the sterile module	I153 225
Battery Li-Ion	for PIPETGIRL	I155 066



Wilford Industrial Estate
Ruddington Lane
Wilford
Nottingham
NG11 7EP
Tel: 0115 982 1111
Email: slsinfo@scientific-labs.com



INTEGRA