

FREE USB Temperature Data Logger when you buy a Lec Refrigerator/Freezer

Lec Medical provide a range of safe, reliable and energy efficient fridges, freezers and fridge-freezers for laboratories. Spark free they are ideal for the safe storage of a wide variety of volatile liquids and chemicals.

All models include:

- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark-free - There is nothing electrical inside the fridge so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces



[See overleaf for ordering information](#)



Laboratory Refrigerators

Code	Model	Capacity, L	Dims, w x d x h, mm	List Price, £
REF4042	LR307C	107	495 x 620 x 830	605.00
REF4044	LSR151	151	595 x 595 x 845	540.00
REF4046	LSR288	288	595 x 634 x 1570	770.00
REF4048	LR1607C	475	695 x 910 x 1688	3,200.00

Laboratory Freezers

Code	Model	Capacity, L	Dims, w x d x h, mm	List Price, £
REF4052	ISU37C	102	495 x 620 x 830	645.00
REF4054	LSF151	151	595 x 595 x 845	595.00
REF4056	LSF232	232	595 x 634 x 1570	820.00

Laboratory Fridge/Freezers

Code	Model	Capacity, L	Dims, w x d x h, mm	List Price, £
REF4062	LSC119	106/13	540 x 600 x 860	425.00
REF4064	LSC263	159/104	540 x 600 x 1680	820.00
REF4066	LSC324	190/134	600 x 600 x 1870	1,170.00

To claim your free USB Temperature Datalogger quote
LOGGER when placing your order