

NEW

Haier -86°C ULT Freezers

- Haier freezers are designed and manufactured to provide safe, reliable, energy efficient long term storage of biological samples.
- Reliable temperature performance -86°C +/-6°C
- Energy Saving VIP insulation and independent inner doors reduce energy consumption without reducing storage capacity
- Microprocessor digital control with 1°C increments, temperature adjustable from -40°C to -86°C

Safety Control System

- Alarms raised if voltage or chamber/ambient temperatures exceed range, or in the event of power failure or sensor error
- Audible and visible alarm indicator
- Multiple protection functions: password, start-up delay and high/low voltage protection

Refrigeration System

- Optimised cascade refrigeration technology
- Patented sealing design keeps the cabinet from frosting
- CO₂ and LN₂ back up systems available

Ergonomic Design

- LED displays chamber, set point and ambient temperatures, plus input voltage, temperature alarm set point and alarm warning.
- Safe door lock to prevent the unauthorized opening.
- Innovative integrated door handle and lock; compacted casters

5 year warranty



Code	Description	Capacity, L	Dims, w x d x h, mm	Shelves / Inner doors	List Price, £
SLS1052	Upright Freezer	338	812 x 893 x 1846	3/2	5,075.00
SLS1054	Upright Freezer	388	812 x 893 x 1980	3/2	5,645.00
SLS1056	Upright Freezer	486	945 x 900 x 1980	3/4	6,325.00
SLS1058	Upright Freezer	490	860 x 900 x 1980	3/4	6,640.00
SLS1060	Upright Freezer	628	1035 x 900 x 1980	3/4	7,085.00
SLS1062	Upright Freezer	728	1145 x 900 x 1980	3/4	7,830.00
SLS1064	Upright Freezer	828	1145 x 980 x 1980	3/4	8,395.00
SLS1066	Chest Freezer	100	769 X 825 x 1120	-/1	4,355.00
SLS1068	Energy Efficient Upright Freezer	578	895 x 980 x 1980	3/4	8,085.00
SLS1070	Energy Efficient Upright Freezer	728	1041 x 980 x 1980	3/4	8,545.00